



Recombinant Arabidopsis thaliana Vesicle-associated protein 4-2 (PVA42)

Product Code	CSB-BP851890DOA
Abbreviation	PVA42
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VYN2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTMTEEKPTS DGKGGWGFFK IPFRNSSGHR NAASSAATSP FPSGASSSST SSHLHNHHQH HHQHHHQHHH QLGYNPHGD GSGQNQHPTP SPSVSSVAKS FLPTKRRLKL DPSEKLYFPY EPGKQVRSI KIKNTSKSHV AFKFQTTAPK SCFMRPPGAI LAPGETIAT VFKFVEPPEN NEKPMQDRSR VKFKIMSLKV KGPMDYVPEL FDEQKDDVSK EQILRVIFLD PERSNPALK LKRQLAEADA AVEARKKPE ETGPKMIGEG LVIDEWKERR ERYLAQQQGE GADSV
Source	Baculovirus
Target Names	PVA42
Protein Names	Recommended name: Vesicle-associated protein 4-2 Alternative name(s): Plant VAP homolog 4-2 Short name= AtPVA42 VAMP-associated protein 4-2
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.