



# Recombinant Human Nucleolar complex protein 3 homolog (NOC3L)

<b>Product Code</b>	CSB-BP851960HU
<b>Abbreviation</b>	NOC3L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WTT2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKARRNKKQI PSFRKLIKTS KVKLENKLN KQFKQQSTLK KYRKEQRKLR QAVKDAVSKK PIPLNPKEK RPGKRIREE EEEEEALPLD MMEDDLQLM KDLGQRVSFL TRDLSSEP V HAKKRKHERI IDKYEKIPRT LQTAPEKELI HLLPIKDKSG IIPQTREKPV TDSNKDEEDQ EEERELEEEI IEDPIQELTI EEHLIERKKK LQEKKMHIAA LASAILSDPE NNIKKLKELR SMLMEQDPDV AVTVRKLIV SLMELFKDIT PSYKIRPLTE AEKSTKTRKE TQKLREFEEG LVSQYKFYLE NLEQMVKDWK QRKLKSNV SLKAYKGLAE VAVKSLCELL VALPHFNFHN NIIVLIVPLM NDMSKLISEM CCEAVKKLFK QDKLGQASLG VIKVISGFVK GRNYEVRPEM LKTFLCLRIK EVEVKKDTEG INKPKKMTF KEKRKSLSRM QRKWKAEEK LERELREAEA SESTEKKLKL HTETLNIVFV TYFRILKKAQ RSPLLPAVLE GLAKFAHLIN VEFFD DLLV LHTLIESGDL SYQESLHCVQ TAFHILSGQG DVLNIDPLKF YTHLYKTLFK LHAGATNEGV EIVLQCLDVM LTKRRKQVSQ QRALAFIKRL CTLALHVLPN SSIGILATTR ILMHTFPKTD LLLDSESQGS GVFLPELDEP EYCNAQNTAL WELHALRRHY HPIVQRFAAH LIAGAPSEGS GALKPELSRR SATELFEAYS MAEMTFNPPV ESSNPKIKGK FLQGDSFLNE DLNQLIKRYS SEVATESPLD FTKYLKTS LH
<b>Source</b>	Baculovirus
<b>Target Names</b>	NOC3L
<b>Protein Names</b>	Recommended name: Nucleolar complex protein 3 homolog Short name= NOC3 protein homolog Alternative name(s): Factor for adipocyte differentiation 24 NOC3-like protein Nucleolar complex-associated protein 3-like protein
<b>Expression Region</b>	1-800
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.