



Recombinant Human Nucleolar complex protein 3 homolog (NOC3L)

Product Code	CSB-YP851960HU
Abbreviation	NOC3L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WTT2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKARRNKKQI PSFRKLIKTS KVKLENKLN KQFKQQSTLK KYRKEQRKLR QAVKDAVSKK PIPLNPKEK RPGKRIEREE EEEEEALPLD MMEDDLQLM KDLGQRVSFL TRDLSSEPVE HAKKRKHERI IDKYEKIPRT LQTAPEKELI HLLPIKDKSG IIPQTREKPV TDSNKDEEDQ EEERELEEEI IEDPIQELTI EEHLIERKKK LQEKKMHIAA LASAILSDPE NNIKKLKELR SMLMEQDPDV AVTVRKLIV SLMELFKDIT PSYKIRPLTE AEKSTKTRKE TQKLREFEEG LVSQYKFYLE NLEQMVKDWK QRKLKSNV SLKAYKGLAE VAVKSLCELL VALPHFNFHN NIIVLIVPLM NDMSKLISEM CCEAVKKLFK QDKLGQASLG VIKVISGFVK GRNYEVRPEM LKTFLCLRIK EVEVKKDTEG INKPKKFTF KEKRKSLSRM QRKWKAEEK LERELREAEA SESTEKKLKL HTETLNIVFV TYFRILKKAQ RSPLLPAVLE GLAKFAHLIN VEFFDILLV LHTLIESGDL SYQESLHCVQ TAFHILSGQG DVLNIDPLKF YTHLYKTLFK LHAGATNEGV EIVLQCLDVM LTKRRKQVSQ QRALAFIKRL CTLALHVLPN SSIGILATTR ILMHTFPKTD LLLDSESQGS GVFLPELDEP EYCNAQNTAL WELHALRRHY HPIVQRFAAH LIAGAPSEGS GALKPELSRR SATELFEAYS MAEMTFNPPV ESSNPKIKGK FLQGDSFLNE DLNQLIKRYS SEVATESPLD FTKYLKTS LH
Source	Yeast
Target Names	NOC3L
Protein Names	Recommended name: Nucleolar complex protein 3 homolog Short name= NOC3 protein homolog Alternative name(s): Factor for adipocyte differentiation 24 NOC3-like protein Nucleolar complex-associated protein 3-like protein
Expression Region	1-800
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.