



Recombinant Human Uncharacterized protein C17orf63 (C17orf63)

Product Code	CSB-BP851963HU
Abbreviation	C17orf63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WU58
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLACLPGPGD LSFQLLSHTQ MNTGLQKWDT TQKMRTAHYP TPAELDAYAK KVANNPLTIK IFPNSVKVPQ RKHVRRTVNG LD TSAQRYSP YPTQAATKAG LLAIVKVPK SILKDFDGTR ARLLEAIMN PPVAPYATVA PSTLAHPQAQ ALARQQALQH AQTLAHAPPQ TLQHPQGIPP PQALSHPQSL QQPQGLGHPQ PMAQTQGLVH PQALAHQGLQ HPHNPLLHGG RKMPDSDAPP NVTVSTSTIP LSMAATLQHS QPPDLSSIVH QINQFCQTRA GISTTSVCEG QIANPSPISR SLLINASTRV STHSVPTPMP SCVVNPMEHT HAATAALPAA GPNLPTGIS RVPTGYPSDL KPVTWNQHQQL AHLQCMCSEA SGTPAPGLTG KHAAGRELAG PGFVGKAPAY PQELCLAQSF HLKPPLEKPT PSPPVNGMAA PLAYPNGHYF QPLWNNILPT PNSDSSGSQD LAMPFHGGQP TGAPLDCAAA PGAHYRAGTG GGPVASQNSL MQTVDYLSGD FQQACFREQS LAMLKAHRA PGNRAPDPTTE SRSLHIQHPG YR
Source	Baculovirus
Target Names	FAM222B
Protein Names	Recommended name: Uncharacterized protein C17orf63
Expression Region	1-562
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.