



# Recombinant Human Uncharacterized protein C17orf63 (C17orf63)

<b>Product Code</b>	CSB-EP851963HU
<b>Abbreviation</b>	C17orf63
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WU58
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLACLPGPGD LSFQLLSHTQ MNTGLQKWDT TQKMRTAHYP TPAELDAYAK KVANNPLTIK IFPNSVKVPQ RKHVRRTVNG LD TSAQRYSP YPTQAATKAG LLAIVKVPK SILKDFDGTR ARLLEAIMN PPVAPYATVA PSTLAHPQAQ ALARQQALQH AQTLAHAPPQ TLQHPQGIPP PQALSHPQSL QQPQGLGHPQ PMAQTQGLVH PQALAHQGLQ HPHNPLLHGG RKMPDSDAPP NVTVSTSTIP LSMAATLQHS QPPDLSSIVH QINQFCQTRA GISTTSVCEG QIANPSPISR SLLINASTRV STHSVPTPMP SCVVNPMEHT HAATAALPAA GPNLPTGIS RVPTGYPSDL KPVTWNQHL AHLQMCSEA SGTPAPGLTG KHAAGRELAG PGFVGKAPAY PQELCLAQSF HLKPPLEKPT PSPPVNGMAA PLAYPNGHYF QPLWNNILPT PNSDSSGSQD LAMPFHGGQP TGAPLDCAAA PGAHYRAGTG GGPVASQNSL MQTVDYLSGD FQQACFREQS LAMLKAHRA PGNRAPDPT ESRSLHIQHPG YR
<b>Source</b>	E.coli
<b>Target Names</b>	FAM222B
<b>Protein Names</b>	Recommended name: Uncharacterized protein C17orf63
<b>Expression Region</b>	1-562
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.