



Recombinant Human Sorting nexin-33 (SNX33)

Product Code	CSB-MP851969HU
Abbreviation	SNX33
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WV41
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MALKGRALYD FHSENKEEIS IQQDEDLVIF SETSLDGWLQ GQNSRGETGL FPASYVEIVR SGISTNHADY SSSPAGSPGA QVSLYNPSV ASPARSGGGS GFLSNQGSFE EDDDDDDWDDW DDGCTVVEEP RAGGLGTNGH PPLNLSYPGA YPSQHMAFRP KPPLERQDSL ASAKRGSVVG RNLNRFSCFV RSGVEAFILG DVPMMAKIAE TYSIEMGPRG PQWKANPHPF ACSVEDPTKQ TKFKGIKSYI SYKLTPTHAA SPVYRRYKHF DWLYNRLLHK FTVISVPHLP EKQATGRFEE DFIEKRKRRL ILWMDHMTSH PVLSQYEGFQ HFLSCLDDKQ WKMGKRRAEK DEMVGASFL TFIQIPTEHQD LQDVEDRVDT FKAFSKKMDD SVLQLSTVAS ELVRKHVGGF RKEFQKLGS FQAISHSFQM DPPFCSEALN SAISHTGR TY EAIGEMFAEQ PKNDLFQMLD TSLYQGLLS NFPDIIHLQK GAFAKV KESQ RMSDEGRMVQ DEADGIRRRRC RVVGFALQAE MNHFHQRREL DFKHMMQNYL RQQLFYQRV GQQLEKTLRM YDNL
Source	Mammalian cell
Target Names	SNX33
Protein Names	Recommended name: Sorting nexin-33 Alternative name(s): SH3 and PX domain-containing protein 3
Expression Region	1-574
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.