



# Recombinant Human E3 ubiquitin-protein ligase RNF138 (RNF138)

<b>Product Code</b>	CSB-YP851971HU
<b>Abbreviation</b>	RNF138
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WVD3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AEDLSAATS YTEDDFYCPV CQEVLKTPVR TTACQHVFCR KCFLTAMRES GAHCPLCRGN VTRRERACPE RALDLENIMR KFSGSCRCCA KQIKFYRMRH HYKSCCKYQD EYGVSSIIPN FQISQDSVGN SNRSETSTSD NTETYQENTS SSGHPTFKCP LCQESNFTRQ RLLDHCNSNH LFQIVPTCP ICVSLPWGDP SQITRNFVSH LNQRHQFDYG EFNVLQLDEE TQYQTAVEES FQVNI
<b>Source</b>	Yeast
<b>Target Names</b>	RNF138
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase RNF138 EC= 6.3.2.- Alternative name(s): Nemo-like kinase-associated RING finger protein Short name= NLK-associated RING finger protein Short name= hNARF RING finger protein 138
<b>Expression Region</b>	2-245
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. Alternatively spliced transcript variants encoding distinct isoforms have been observed.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.