



Recombinant Human Ubiquitin-like domain-containing CTD phosphatase 1 (UBLCP1)

Product Code	CSB-EP851974HU
Abbreviation	UBLCP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WVY7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ALPIIVKWG GQEYSVTTLS EDDTVLDLKQ FLKTLTGVLV ERQKLLGLKV KGKPAENDVK LGALKLKPNT KIMMMGTREE SLEDVLGPPP DNDDVVNDFD IEDEVVEVEN REENLLKISR RVKEYKVEIL NPPREGKKLL VLDVDYTLFD HRSCAETGVE LMRPYLHEFL TSAYEDYDIV IWSATNMKWI EAKMKELGVS TNANYKITFM LDSAAMITVH TPRRGLIDVK PLGVIWGKFS EFYSKNTIM FDDIGRNFLM NPQNGLKIRP FMKAHLNRDK DKELLKLTQY LKEIAKLDDF LDLNHKYWER YLSKKQGQ
Source	E.coli
Target Names	UBLCP1
Protein Names	Recommended name: Ubiquitin-like domain-containing CTD phosphatase 1 EC= 3.1.3.16 Alternative name(s): Nuclear proteasome inhibitor UBLCP1
Expression Region	2-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.