



Recombinant L-threonine 3-dehydrogenase (tdh)

Product Code	CSB-EP852073EOD
Abbreviation	tdh
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8XEJ1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MKALSCLKAE EGIWMTDVPV PELGHNDLLI KIRKTAICGT DVHIYNWDEW SQKTIPVPMV VGHEYVGEVV GIGQEVKGFK IGDRVSGEGH ITCGHCRNCR GGRTHLCRNT IGVGVNRPGC FAEYLVIPAF NAFKIPDNIS DDLASIFDPF GNAVHTALSF DLVGEDVLVS GAGPIGIMAA AVAKHVGARN VVITDVNEYR LELARKMGIT RAVNVAKENL NDVMAELGMT EGFDVGLEMS GAPPAFRTML DTMNHGGRIA MLGIPPSDMS IDWTKVIFKG LFIKGIYGRE MFETWYKMAA LIQSGLDLSP IITHRFSIDD FQKGFDMCS GQSGKVILSW D
Source	E.coli
Target Names	tdh
Protein Names	Recommended name: L-threonine 3-dehydrogenase EC= 1.1.1.103
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>