



Recombinant Chicken Transcription factor E2F1 (E2F1)

Product Code	CSB-EP852673CH-B
Abbreviation	E2F1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q90977
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MATAGGAAGL AALLGGASPH LLIVSASEEP AGGCRPDADL LLFATPQPSR PGPAPRRPAL GRPPVKRKLN LETDHQYIAE SLPAARGRAR IPGRGAKSPG EKSRYETSLN LTTKRFLLELL SQSPDGVVDL NWAAEVLKVQ KRRIYDITNV LEGIQLITKK SKNNIQWLGS QVAAGASSRQ RLLKELRDLD QAAERQLDDL IQTCTVRLRL LTEDPSNQHA AYTVCQDLRS IVDPSEQMVM VIKAPPETQL QVSDPGEAFQ VSVRSTQGPI DVFLCPEDSS GVCSPVKSPF KAPAEELSPG SSQQRASPLL HSAQDVNMLL PEALLPGTAL PTKCPTEDVS LSPLASMDTL LEHGKDDDFPG FLADEFIALS PPQPQDYHFG LEEGEGISEL FDCDFGDFTH LDF
Source	E.coli
Target Names	E2F1
Protein Names	Recommended name: Transcription factor E2F1 Short name= E2F-1
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the E2F family of transcription factors. The E2F family plays a crucial role in the control of cell cycle and action of tumor suppressor proteins and is also a target of the transforming proteins of small DNA tumor viruses. The E2F proteins contain several evolutionally conserved domains found in most members of the family. These domains include a DNA binding domain, a dimerization domain which determines interaction with the differentiation regulated transcription factor proteins (DP), a transactivation domain enriched in acidic amino acids, and a tumor suppressor protein association domain which is embedded within the transactivation domain. This protein and another 2 members, E2F2 and E2F3, have an additional cyclin



binding domain. This protein binds preferentially to retinoblastoma protein pRB in a cell-cycle dependent manner. It can mediate both cell proliferation and p53-dependent/independent apoptosis.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.