



Recombinant Bothrops jararacussu Venom nerve growth factor (NGF)

Product Code	CSB-EP852675BUX-B
Abbreviation	NGF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q90W38
Product Type	Recombinant Protein
Immunogen Species	Bothrops jararacussu (Jararacussu)
Purity	≥85% (SDS-PAGE)
Sequence	EDHPVHNR GEYSVCDSVN VVWANKTTAT DIRGNVVTVM VDVNINNNVY KQYFFETKCR NPNPVPTGCR GIDARHWNSY CTTTNTFVKA LTMEGNQASW RFIRIDTACV CVISRKNENF G
Source	E.coli
Target Names	NGF
Protein Names	Recommended name: Venom nerve growth factor Short name= v-NGF Short name= vNGF Alternative name(s): Bj-NGF
Expression Region	123-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.