



Recombinant *Cnemidophorus uniparens* Estrogen receptor (ESR1)

Product Code	CSB-MP852710CLAZ
Abbreviation	ESR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91424
Product Type	Recombinant Protein
Immunogen Species	<i>Aspidoscelis uniparens</i> (Desert grassland whiptail lizard) (<i>Cnemidophorus uniparens</i>)
Purity	>85% (SDS-PAGE)
Sequence	GHNDYMCPAT NQCTIDKNRR KSCQACRLRK CYEVGMMKGG IRKDRRGGRILKHKRQREEH DNRNAGAIVE RRSPNLWPSP LMITHNKKNS PALSLTADQIVSALLEAEPP VVYSEYDPSR PFSEASMMTL LTNLADRELV HMINWAKRVP GFVDLSLHDQ VHLLCAWLE ILMIGLVWRS VEHPGKLLFA PNLLLDNRNQG KCVEGFVEIF DMLLATSSRF RMMNVQGEEF VCLKSILLN SGIYTFLSST LKSLEEKDHI HRVLDKIIDT LLHLMAKSGL SLQQQHRLA QLLLILSHFR HMSNKGM
Source	Mammalian cell
Target Names	ESR1
Protein Names	Recommended name: Estrogen receptor Short name= ER Alternative name(s): ER-alpha Estradiol receptor Nuclear receptor subfamily 3 group A member 1
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.