



Recombinant Mouse Epsin-3 (Epn3)

Product Code	CSB-EP852779MO-B
Abbreviation	Epn3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91W69
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MTTSALRRQV KNIVHNYSEA EIKVREATSN DPWGPPSSLM SEIADLTFNT VAFAEVMGMV WRRLNDSGKN WRHVYKALT LDYLLKTGSE RVAHQREN YTIQTLKDFQ YIDRDGKDQG VNVREKVKQV MALLKDEERL RQERTHAKT KERMALEGMG IGSGQLGYSR RSRGSPSSYT SASSSPRYAS DLEQARPQTS GEEELQLQLA LAMSREEAER PVPPASHRDE DLQLQLALS SRQEHEKGVR SWKGDDSPVA NGAEPAGQRR QRDREPEREE RKEEEKLKTS QSSILDLADI FAPAPALPST HCSADPWDIP GLRPNTEPSG SSWGPSADPW SPVPSGNALS RSQPWDLPT LSSSEPWGRT PVLPSGPIA DPWAPSSPTR KLPSTGADPW GASMETS DTS ALGGASPFDP FAKPLESTEP KESRDSAQAL PTGKSPSTVE LDPFGDSSPS CKQNGMKEPE ALDLGVLGEA LPQQPGKEAR PCRTPEFLG PSASSLVNLD SLVKAPLAAR TRNPFLTGLG VPSPTNPFGA GDQGRPTLNQ MRTGSPALGL PPGGPVGAPV GSMTYSASLP LPLSSVPVGA TLPASVSVFP QAGAFAPPPA SLPQPLLPTS GPMGPLPPQA GTNPFL
Source	E.coli
Target Names	Epn3
Protein Names	Recommended name: Epsin-3 Alternative name(s): EPS-15-interacting protein 3
Expression Region	1-636
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.