



# Recombinant Mouse Epsin-3 (Epn3)

<b>Product Code</b>	CSB-MP852779MO
<b>Abbreviation</b>	Epn3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q91W69
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTSALRRQV KNIVHNYSEA EIKVREATSN DPWGPPSSLM SEIADLTFNT VAFAEVMGMV WRRLNDSGKN WRHVYKALT LDYLLKTGSE RVAHQCRENL YTIQTLKDFQ YIDRDGKDQG VNVREKVKQV MALLKDEERL RQERTHALKT KERMALEGMG IGSGQLGYSR RSRGSPSSYT SASSSPRYAS DLEQARPQTS GEEELQLQLA LAMSREEAER PVPPASHRDE DLQLQLALSL SRQEHEKGVR SWKGDSPVA NGAEPAGQRR QRDREPEREE RKEEEKLKTS QSSILDLADI FAPAPALPST HCSADPWDIP GLRPNTEPSG SSWGPSADPW SPVPSGNALS RSQPWDLPT LSSSEPWGRT PVLPSGPIA DPWAPSSPTR KLPSTGADPW GASMETSPTS ALGGASPFDP FAKPLESTEP KESRDSAQAL PTGKSPSTVE LDPFGDSSPS CKQNGMKEPE ALDLGVLGEA LPQQPGKEAR PCRTPEPFLG PSASSLVNLD SLVKAPLAAR TRNPFLTGLG VPSPTNPFGA GDQGRPTLNQ MRTGSPALGL PPGGPVGAPV GSMTYSASLP LPLSSVPVGA TLPASVSVFP QAGAFAPPA SLPQPLLPTS GPMGPLPPQA GTNPFL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Epn3
<b>Protein Names</b>	Recommended name: Epsin-3 Alternative name(s): EPS-15-interacting protein 3
<b>Expression Region</b>	1-636
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.