



Recombinant Mouse Diacylglycerol kinase gamma (Dgkg)

Product Code	CSB-YP852784MO
Abbreviation	Dgkg
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91WG7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSEEQWVSLSEEFDQLQKYSEYSSKKIKDVLAEFNEGGS LRQYDPHKPI SYDVFKLFMR AYLEVDLPQP LSTHLFLAFS QKPRQETPDH PKEGASSSEP NVSDYNSDNA AKADEACAPD TESKTTKTQA PSKELEAAP WEDPGALASS SDAPVVYLKD VVCYLSLMET GRPQDKLEFM FRLYDSDENG LLDQAEMDQI VSQMLHVAQY LEWDPTELRP ILKEMLQGM YDKDGFVSLQ EWINGGMTTI PLLVLLGMDD SGSKGDGRHA WTLKHFKKPT YCNFCRAMLM GVGKQGLCCI YCKYTVHQRC VSKTIHGCVK TNSKAKRSGE VMQHAWVEGN SSVKCDRCHK SIKCYQSVTA RHCVWCRMTF HRKCELSTVC DGGELKDHIL LPTSICPVSG DRQGGKSDGS VAAKGELVTQ YKIIPSPGTH PLLVLVNPKS GGRQGERILR KFHLLNPEQ VFNLNNGGPT PGLNFFHDTP DFRVLACGGD GTVGWILDCI DKANFTKHPP VAVLPLGTGN DLARCLRWGG GYEGGSLTKI LKEIEQSPLV MLDRWYLEVM PREEVENDQ VPYNIMNNYF SIGVDASIAH RFHMMREKHP EKFNARMKNK LWYFEFGTSE TFAATCKKLH DHIELECDGV EVDLSNIFLE GAILNIPSM YGGTNLWGET KKNRAVIRE RKSVDPKEL KCCVQDLSAQ LLEVVGLEGA MEMGQIYTGL KSAGRRLAQC SSVTIRTNKL LPMQVDGEPW MQPQCTIKIT HKNQAPMMM PPQKSSFFSL RRKSRSKD
Source	Yeast
Target Names	Dgkg
Protein Names	Recommended name: Diacylglycerol kinase gamma Short name= DAG kinase gamma EC= 2.7.1.107 Alternative name(s): 88 kDa diacylglycerol kinase Diglyceride kinase gamma Short name= DGK-gamma
Expression Region	1-788
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes an enzyme that is a member of the type I subfamily of



diacylglycerol kinases, which are involved in lipid metabolism. These enzymes generate phosphatidic acid by catalyzing the phosphorylation of diacylglycerol, a fundamental lipid second messenger that activates numerous proteins, including protein kinase C isoforms, Ras guanyl nucleotide-releasing proteins and some transient receptor potential channels. Diacylglycerol kinase gamma has been implicated in cell cycle regulation and in the negative regulation of macrophage differentiation in leukemia cells. Multiple transcript variants encoding different isoforms have been found for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.