



Recombinant Mouse RING finger protein 37 (Ubox5)

Product Code	CSB-MP852887MO
Abbreviation	Ubox5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q925F4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MVVNLCLPQF RPRIHCNKVS ADGYEVENLI SEDLIKRSYG FRTEYFIRPP IYVTVSFPFN VEICRVNIDL TTGGYQNVSG LELYTSALSS RVSQDAQDCW TTGPVETSVP DKEAFTLVGK VLLKNQNHVV FSHRGFKARP PFSPMEVTLL SPAVVAQELW NKGALSLSHV AHLKIGITHV TGSGISCIKR LEVWGQPART CSQEVINSVL LIASESLPQD LDHAPALPM ESDCDPGGQS ESQHSPCTLQ DMSEVESDVP EEFLDPITLE IMPCPMLLPs GKVIDQSTLE KCNLSEAAWG RVSPDPFTGL AFTPQSQPLP HPSLKARIDR FLLQHSISGC RLLGRAQTPS AMTPSVITLP SRKRKTEQAE HSSHYSLGMS ASSSATSPLF SPTTSEPTAK KMKATSELGL TDMDCSAGPV SHEQKLAQSL EIALTSTLGS MPSFTARLTK GQLQLGTRGS SACRRPASSS EHPRSVSGPE CASCKQAFSS YSTNEPVYQL PCGHLLCRPC LSEKQRSQPM MCTACRQPVT SQDVLRVHF
Source	Mammalian cell
Target Names	Ubox5
Protein Names	Recommended name: RING finger protein 37 Alternative name(s): U-box domain-containing protein 5 UbcM4-interacting protein 5 Ubiquitin-conjugating enzyme 7-interacting protein 5
Expression Region	1-539
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.