



Recombinant Human Protein TFG (TFG)

Product Code	CSB-BP852902HU
Abbreviation	TFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92734
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNGQLDLSGK LIIKAQLGED IRRIPHNED ITYDELVLMM QRVFRGKLLS NDEVTIKYKD EDGDLITIFD SSDLSFAIQC SRILKLTFLV NGQPRPLESS QVKYLRRELI ELRNKVNRL DSLEPPGEPG PSTNIPENDT VDGREEKSAS DSSGKQSTQV MAASMSAFDP LKNQDEINKN VMSAFGLTDD QVSGPPSAPA EDRSGTPDSI ASSSSAAHPP GVQPQQPPYT GAQTQAGQIE GQMYQQYQQQ AGYGAQQPQA PPQQPQQYGI QYSASYSQQT GPQQPQQFQG YGQQPTSQAP APAFSGQPQQ LPAQPPQQYQ ASNYPQTYT AQTSQPTNYT VAPASQPGMA PSQPGAYQPR PGFTSLPGST MTPPPSGPNP YARNRPPFGQ GYTQPGPGYR
Source	Baculovirus
Target Names	TFG
Protein Names	Recommended name: Protein TFG Alternative name(s): TRK-fused gene protein
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.