



Recombinant Rickettsia conorii Putative uncharacterized protein RC0841 (RC0841)

Product Code	CSB-MP853000RMS
Abbreviation	RC0841
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92HD0
Product Type	Recombinant Protein
Immunogen Species	Rickettsia conorii (strain ATCC VR-613 / Malish 7)
Purity	>85% (SDS-PAGE)
Sequence	MQLLLLIPLH YRLKDSKRIY QIILNLGSLI RISSIRQDNM ELLELIIDPM IKAEVIERVS EFLTKTLGKI IVKCNDSGFI IANRVGYFLL ELVARKAISQ NLDVATIDKI FTTFLGLPST GIFGLYDLIG YDVMKLISSS LLASLPANDA YHKIICKDPC S
Source	Mammalian cell
Target Names	RC0841
Protein Names	Recommended name: Putative uncharacterized protein RC0841
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.