



Recombinant *Oryza sativa* subsp. *japonica* Protein disulfide isomerase-like 2-2 (PDIL2-2)

Product Code	CSB-EP853186OFG
Abbreviation	PDIL2-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q942L2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	APA FADGDDVVAL TESTFEKEVG QDRGALVEFY APWCGHCKKL APEYEKLGAS FKKAKSVFIA KVDCDEHKS SV CSKYGVSGYP TIQWFPKGSL EPKKYEGQRS AEALAEFVNT EGGTNVKLAT IPSSVVVLGP DNFDIVLDE NKDILVEFYA PWCGHCKHLA PIYEKLSVY KLDDGVVIAN LDADKHKDLA EKYGVSGYPT LKFFPKGNKA GEDYDGGREL DDFVKFINEK CGTSRDTKGG LTSEAGRIAS LDALAKEFLG AANDKRKEIL SNMEEEVVKL SGSAAKHGKV YIAIAKKILD KGHDYTKKET ERLERMLEKS ISPSKADEFI IKKNVLSTFS S
Source	E.coli
Target Names	PDIL2-2
Protein Names	Recommended name: Protein disulfide isomerase-like 2-2 Short name= OsPDIL2-2 EC= 5.3.4.1 Alternative name(s): Protein disulfide isomerase-like 4-2 Short name= OsPDIL4-2
Expression Region	28-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.