



Recombinant Arabidopsis thaliana

Uncharacterized protein At2g17340 (At2g17340)

Product Code	CSB-EP853208DOA-B
Abbreviation	At2g17340
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q949P3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MESDSEMVPF PQLPMPIENN YRACTIPYRF PSDDPKKATP NEISWINVFA NSIPSFKKRA ESDITVPDAP ARAEKFAERY AGILEDLKKD PESHGGPPDG ILLCRLREQV LRELGFRDIF KVKKDEENAK AISLFPQVVS LSDAIEDDGK RLENLVRGIF AGNIFDLGSA QLAEVFSRDG MSFLASCQNL VPRPWVIDDL ENFQAKWINK SWKKAVIFVD NSGADIILGI LPFARELLRR GAQVVLAANE LPSINDITCT ELTEILSQLK DENGQLLQVD TSKLLIANSG NDLPVIDLSR VSQELAYLSS DADLVIVEGM GRGIETNLYA QFKCDSLKIG MVKHLEVAEF LGGRLYDCVF KFNEVQS
Source	E.coli
Target Names	At2g17340
Protein Names	Recommended name: Uncharacterized protein At2g17340
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.