



Recombinant Human General transcription factor 3C polypeptide 6 (GTF3C6)

Product Code	CSB-EP853368HU-B
Abbreviation	GTF3C6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q969F1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAAADERSP EDGEDEEEEE QLVLVELSGI IDSDFLSKCE NKCKVLGIDT ERPILQVDSC VFAGEYEDTL GTCVIFEENV EHADTEGNNK TVLKYKCHTM KKLSMTRLL TEKKEGEENI GGVEWLQIKD NDFSYPNMI CNFLHENEDE EVVASAPDKS LELEEEIIM NDSSNLSCEQ EKPMHLEIED SGPLIDIPSE TEGSVFMETQ MLP
Source	E.coli
Target Names	GTF3C6
Protein Names	Recommended name: General transcription factor 3C polypeptide 6 Alternative name(s): Transcription factor IIIC 35 kDa subunit Short name= TFIIIC 35 kDa subunit Short name= TFIIIC35 Transcription factor IIIC subunit 6
Expression Region	2-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.