



Recombinant Human Leukocyte receptor cluster member 9 (LENG9)

Product Code	CSB-MP853388HU
Abbreviation	LENG9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96B70
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGSRPPCGAT SSARRACQFP APMAAAREPE LPQEAPATEP APPPACRFFL EGRCRFGARC RQPHPGAPAP PGREAQPEAG AKKPLRTAA DVIQRIRWDP RLDPADFSVG YVDRFLGVRE EPFSAFCWDQ PLAALGPGVL AVPQHRVRF RFHGRLVWDR ASRTDLVFGS GSAAGRGPTI LDAPNTEGAH GAEGAEWTLA GTGQEAQAAP KRGSTRPLCT GHQEPGVEEP GELEAAQERA LGTAADLGT APRGRLAGVT EEALKPTAAT RTLLGGKEA QALGVPGGSA ETTEAEWGPA AWPEDKRARL SVAAPCQPRP THFVALMVTE PGLQAEVTKA QEYLVHVAPH CANFLVPSQN LHLTLALLRL AGAGEEAAAI GALRRALLAP GLNAPPRLSF RKLVLLGPHV LCAPPSPTLE SMAQVLSQRL EA EGLSTLQS PGQLHPHLTV AKVPHGSQVH LPKLEFTLSQ EVGCQPLQTL WLCRIGRTGG PFQPLAEIRL E
Source	Mammalian cell
Target Names	LENG9
Protein Names	Recommended name: Leukocyte receptor cluster member 9
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.