



Recombinant Human Protein FAM105B (FAM105B)

Product Code	CSB-EP853391HU-B
Abbreviation	FAM105B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BN8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSRGTMPQPE AWPGASCAET PAREAAATAR DGGKAAASGQ PRPEMQCPAE HEEDMYRAAD EIEKEKELLI HERGASEPRL SVAPEMDIMD YCKKEWRGNT QKATCMKMGY EEVSQKFTSI RRVRGDNYCA LRATLFQAMS QAVGLPPWLQ DPELMLLPEK LISKYNWIKQ WKLGLKFDGK NEDLVDKIKE SLTLLRKKWA GLAEMRTAEA RQIACDELFT NEAEEYSLYE AVKFLMLNRA IELYNDKEKG KEVPPFSVLL FARDTSNDPG QLLRNHLNQV GHTGGLEQVE MFLLAYAVRH TIQVYRLSKY NTEEFITVYP TDPPKDWPVV TLIAEDDRHY NIPVRVCEET SL
Source	E.coli
Target Names	OTULIN
Protein Names	Recommended name: Protein FAM105B
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.