



Recombinant Human Cob (I)yrinic acid a,c-diamide adenosyltransferase, mitochondrial (MMAB)

Product Code	CSB-MP853413HU
Abbreviation	MMAB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96EY8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QSRGPQGV EDGDRPQPSS KTRIPRIKIYT KTGDKGFSST FTGERRPKDD QVFEAVGTTD ELSSAIGFAL ELVTEKGHTF AEELQKIQCT LQDVGSALAT PCSSAREAHL KYTTFKAGPI LELEQWIDKY TSQLPPLTAF ILPSGGKISS ALHFCRAVCR RAERRVVPLV QMGETDANVA KFLNRLSDYL FTLARYAAMK EGNQEKIYMK NDPSAESEGL
Source	Mammalian cell
Target Names	MMAB
Protein Names	Recommended name: Cob(I)yrinic acid a,c-diamide adenosyltransferase, mitochondrial EC= 2.5.1.17 Alternative name(s): Cob(I)alamin adenosyltransferase Methylmalonic aciduria type B protein
Expression Region	33-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a protein that catalyzes the final step in the conversion of vitamin B(12) into adenosylcobalamin (AdoCbl), a vitamin B12-containing coenzyme for methylmalonyl-CoA mutase. Mutations in the gene are the cause of vitamin B12-dependent methylmalonic aciduria linked to the cblB complementation group.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.