



Recombinant Human AMSH-like protease (STAMBPL1)

Product Code	CSB-EP853417HU-B
Abbreviation	STAMBPL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FJ0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDQPFTVNSL KKLAAMPDHT DVLSPEERV RALSKLGCNI TISEDITPRR YFRSGVEMER MASVYLEEGN LENAFLVLYNK FITLFVEKLP NHRDYQQCAV PEKQDIMKKL KEIAFPRTDE LKNDLLKKYN VEYQEYLQSK NKYKAEILKK LEHQRLIEAE RKRIAQMRQQ QLESEQFLFF EDQLKKQELA RGQMRSQQTS GLSEQIDGSA LSCFSTHQNN SLLNVFADQP NKSDATNYAS HSPPVNRALT PAATLSAVQN LVVEGLRCVV LPEDLCHKFL QLAESNTVRG IETCGILCGK LTHNEFTITH VIVPKQSAGP DYCDMENVEE LFNVDQHDH LTLGWIHHP TQTAFLSSVD LHTHCSYQLM LPEAIAIVCS PKHKDTGIFR LTNAGMLEVS ACKKKGFHHP TKEPRLFSIC KHVLVKDIKI IVLCLR
Source	E.coli
Target Names	STAMBPL1
Protein Names	Recommended name: AMSH-like protease Short name= AMSH-LP EC= 3.4.19.- Alternative name(s): STAM-binding protein-like 1
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.