



Recombinant Human Protein LTV1 homolog (LTV1)

Product Code	CSB-EP853422HU-B
Abbreviation	LTV1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96GA3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPHRKKKPF I EKKKAVSFHL VHR SQRDPLA ADESAPQRVL LPTQKIDNEE RRAEQRK YGV FFDDDYDYLQ HLKEPSGPSE LIPSS TFSAH NRREEKEETL VIPSTGIKLP SSVFASEFEE DVGLLNKAAP VSGPRLDFDP DIVAALDDDF DFDDPDN LLE DDFILQANKA TGE EEGMDIQ KSENE DDSEW EDVDDEK GDS NDDYDSAGLL SDEDCMSVPG KTHRAIADHL FWSEETKSRF TEYSMTSSVM RRNEQLTLHD ERFEKFYEQY DDDEIGALDN AELEGSIQVD SNRLQEV LND YYKEKAENC V KLNTLEPLED QDLPMNELDE S EEEEMITVV LEEAKEKWDC ESICSTYSNL YNHPQLIKYQ PKPKQIRISS KTG IPLNVLP KKGLTAKQTE RIQMINGS DL PKVSTQPRSK NESKEDKRAR KQAIKEERKE RRVEKKANKL AFKLEKRRQE KELLNLKKNV EGLKL
Source	E.coli
Target Names	LTV1
Protein Names	Recommended name: Protein LTV1 homolog
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.