



# Recombinant Human Protein FAM71F1 (FAM71F1)

<b>Product Code</b>	CSB-MP853446HU
<b>Abbreviation</b>	FAM71F1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96KD3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLSSFPHRKT WRKSKKTVKV TRSYPTFPSL NAWEEFRGLL PVDGEPNPGA GLGVEEGLLC RVVHSPEFNL FLDSVVFESN FIQVKRGRNW RDVYKASNTM ALGVTSSVPC LPLPNILLMA SVKWHQGQONQ TWRNPSIAPN IFLKRILPLR FVELQVCDHY QRILQLRTVT EKIIYYLKLHP DHPETVFHFW IRLVQILQKG LSITTKDPRI LVTHCLVPKN CSSPSGDSKL VQKKLQASQP SESLIQLMTK GESEALSQIF ADLHQQNQLS FRSSRKVETN KNSSGKDSSR EDSIPCTCDL RWRASFTYGE WERENPSGLQ PLSLLSTLAA STGPQLAPPI GNSI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FAM71F1
<b>Protein Names</b>	Recommended name: Protein FAM71F1 Alternative name(s): Protein FAM137A Testis development protein NYD-SP18
<b>Expression Region</b>	1-344
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.