



# Recombinant Sheep Antithrombin-III (SERPINC1)

<b>Product Code</b>	CSB-MP021079SH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P32262
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Ovis aries (Sheep)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	HRSPVEDI CTAKPRDIPV NPMCIYRSPE KKATEGEGSE QKIPGATNRR VWELSKANSH FATAFYQHLA DSKNNNDNIF LSPLSISTAF AMTKLGACNN TLKQLMEVFK FDTISEKTSQ QIHFFFAKLN CRLYRKANKS SELVSANRLF GDKSITFNET YQDISEVVYG AKLQPLDFKG NAEQSRLTIN QWISNKTEGR ITDVIPPQAI DEFTVLVLVN TIYFKGLWKS KFSPENTKKE LFYKADGESC SVPMMYQEGK FRYRRVAEGT QVLELPFKGD DITMVLILPK LEKPLAKVER ELTPDMLQEW LDELTTLLV VHMPHFRIED SFSVKEQLQD MGLEDLFSPE KSRLPGIVAE GRNDLYVSDA FHKAFLEVNE EGSEAAASTV ISIAGRSLNL NRVTFQANRP FLVLIREVAL NTIIFMGRVA NPCVN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SERPINC1
<b>Protein Names</b>	Recommended name: Antithrombin-III Short name= ATIII Alternative name(s): Serp C1
<b>Expression Region</b>	33-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a plasma protease inhibitor and a member of the serpin superfamily. This protein inhibits thrombin as well as other activated serine proteases of the coagulation system, and it regulates the blood coagulation cascade. The protein includes two functional domains: the heparin binding- domain at the N-terminus of the mature protein, and the reactive site domain at the C-terminus. The inhibitory activity is enhanced by the presence of heparin. More than 120 mutations have been identified for this gene, many of which are known to cause antithrombin-III deficiency.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.