



Recombinant Human Uncharacterized protein C8orf48 (C8orf48)

Product Code	CSB-BP853456HU
Abbreviation	C8orf48
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96LL4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAICPELAQT DKSALANLSD ETETLKNSTD EVQTSSSFSS SGGRQSSPLT SGSKLEREKQ TPSLEQGDTQ SELLDYKNYE KKLSKKWINY LKLKDSNFER HQPDTKLPT EITRVSDEELN ALQSYCTMKI NLIHRRGDSK KKTSSRHKKL HLGLDVEASE RDAFSCVTPD ELLNRIYFKN MRTTPKQEA AKQHISYQCP YCNRKRAELA LSAFLKQKKT LLESFLLQER IDEHLHTKDF LTRIGEAHQD FPRLSDDPRI IWKRLTEKSH IRYSGFERSE TEQKLQRDGN SACHLPFSLP FLKRLTLIKP ELVIVNDNV
Source	Baculovirus
Target Names	C8orf48
Protein Names	Recommended name: Uncharacterized protein C8orf48
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.