



Recombinant Human Melanoma-associated antigen B10 (MAGEB10)

Product Code	CSB-YP853459HU
Abbreviation	MAGEB10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96LZ2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPRGQKSKLR AREKRRQARG GLEDLIDALD ILEEEEEESPP SASACLKDVF QSSLDGASNN PHGLREAQST STSATAASHT RHPEGVNDQM EERPICTQDL EATDSFPRGP VDEKVILVH YLLYKYQMKE PITKADMLRN VTQMSKSQFP VILSRASEHL ELIFGLDLKE VEPNKHIYVL VNKLDLGCD A KLSDETGVPK TGLLMTVLGI IFTNGNCVAE EEVWKVFNTM GLYDGIIEHFM FGEPKLLTK DLVKENYLEY QQVPNSDPPR YQFLWGPRAH AETSKMKVLE FLAKVNDTAP SEFSNWyTEA LQDEEERARA RVAAKARVSA TAGARSKVKS SKSSQLQ
Source	Yeast
Target Names	MAGEB10
Protein Names	Recommended name: Melanoma-associated antigen B10 Alternative name(s): MAGE-B10 antigen
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.