



Recombinant Human Melanoma-associated antigen B18 (MAGEB18)

Product Code	CSB-BP853461HU
Abbreviation	MAGEB18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M61
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPRGQKSKLR AREKRHQARC ENQDLGATQA TVAEGESPPS AYLLFGDRPQ NLPAAETPSI PEALQGAPST TNAIAPVSCS SNEGASSQDE KSLGSSREAE GWKEDPLNKK VVSLVHFLQ KYETKEPITK GDMIKFVIRK DKCHFNEILK RASEHMELAL GVDLKEVDPI RHYYAFFSKL DLTYDETTSD EEKIPKTGLL MIALGVIFLN GNRAPPEAVW EIMNMMGVYA DRKHFLYGD RKVMTKDLVQ LKYLEYQQVP NSDPPRYEFL WGPRAHAETS KMKVLEFVAK IHDTVPSAFP SCYEEALRDE EQRTQARAAA RAHTAAMANA RSRTTSSSFS HAK
Source	Baculovirus
Target Names	MAGEB18
Protein Names	Recommended name: Melanoma-associated antigen B18 Alternative name(s): MAGE-B18 antigen
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.