



# Recombinant Human Zinc finger SWIM domain-containing protein 3 (ZSWIM3)

<b>Product Code</b>	CSB-BP853466HU
<b>Abbreviation</b>	ZSWIM3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96MP5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MELGSCFKTY EDFKECF SAY KRENRC SFIL RDCVSVRFHN LNHGTSIRED ILYVQVKFVC IRTQSNRKRT READMCPAYL LRYNERLDR LFISELNTQH IHGDSKVASP GGD TTGKSQK TMCLQRLQPV QPTTKKDLDT AEKSLVEPSF CLDKVQVSSK PEQEGITPSD LAKIAKVMKN FLKVDEGSMA SFSVGDSQHL DRLSFQSSKM TDLFIRFPEN LLLHRVENTQ GHILYAFLVE NKERESRVVH FAVLKAETVT SVAKMLSIFT EFNSDWPKVK VVFVDP SFHY RAILQEIFPA ARILLSIYHT TRLLEK LHR SSANPSFKRL MKEALREAVF VTSEASLKNL CQMSQAVLDE DLFNFLQAHW FTCELLWYMH VRKGLLACNT YMDSL DIVTS KVSSLFREQQ SLLDCILCFV DYIDFFNTKG LKNLPTPPPK LKRARPASMP LKSKKAFGIC GESLTSLPAE ETKPDAQVQ VQQQSQVPPS QVGMLD TLHQ SGSELAYKLC HNEWVVQNS THLVDMAGSS VDVQLLED SH QVSKDGCSCS CSFQQWYHLP CRHILALLHT SQQPVG EAMV CRRWQKKYQY LLGPN GELQD RGMVPNTGQP EKQGRNDMIQ DLSRELANLL MQTEGPELEE RYSTLRKIVD IWAGPSQPSE LFQQPGDFKD VGRLPFLWGK QEEGEGFPPA TAVMHY
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZSWIM3
<b>Protein Names</b>	Recommended name: Zinc finger SWIM domain-containing protein 3
<b>Expression Region</b>	1-696
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.