



Recombinant Human Apoptosis-inducing factor 3 (AIFM3)

Product Code	CSB-MP853472HU
Abbreviation	AIFM3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NN9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGGCFSKPKP VELKIEVVLP EKERGKEELS ASGKGSPRAY QGNGTARHFH TEERLSTPHP YPSPQDCVEA AVCHVKDLEN GQMREVELGW GKVLLVKDNG EFHALGHKCP HYGAPLVKGV LSRGRVRCPW HGACFNISTG DLEDFPGLDS LHKFQVKIEK EKVYVRASKQ ALQLQRRTKV MAKCISPSAG YSSSTNVLIV GAGAAGLVCA ETLRQEGFSD RIVLCTLDRH LPYDRPKLSK SLDTQPEQLA LRPKEFFRAY GIEVLTEAQQV VTDVVRTKKV VFKDGFKLEY SKLLLAPGSS PKTLSCKGKE VENVFTIRTP EDANRVVRLA RGRNVVVVGA GFLGMEVAAY LTEKAHSVSV VELEETPFRR FLGERVGRAL MKMFENNRVK FYMQTEVSEL RGQEGKLKEV VLKSSKVVRA DVCVVGIGAV PATGFLRQSG IGLDSRGFIP VNKMMQTNVP GVFAAGDAVT FPLAWRNNRK VNIPHWQMAH AQGRVAAQNM LAQEAEMSTV PYLWTAMFGK SLRYAGYGEG FDDVIIQGDL EELKFVAFYT KGDEVIAVAS MNYDPIVSKV AEVLASGRAI RKREVELFVL HSKTGDMSWL TGKGS
Source	Mammalian cell
Target Names	AIFM3
Protein Names	Recommended name: Apoptosis-inducing factor 3 EC= 1.-.-.- Alternative name(s): Apoptosis-inducing factor-like protein
Expression Region	1-605
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.