



Recombinant Human Coiled-coil domain-containing protein 115 (CCDC115)

Product Code	CSB-YP853473HU
Abbreviation	CCDC115
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NT0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAALDLRAEL DSLVLQLLGD LEELEGKRTV LNARVEEGWL SLAKARYAMG AKSVGPLQYA SHMEPQVCLH ASEAQEGLQK FKVVRAGVHA PEEVGPREG LRRRKGPSTK PEPESSEAPQ DPLNWFILV PHSLRQAQAS FRDGLQLAAD IASLQNRIDW GRSQLRGLQE KLIKLEPGAA
Source	Yeast
Target Names	CCDC115
Protein Names	Recommended name: Coiled-coil domain-containing protein 115
Expression Region	1-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.