



# Recombinant Human Protein FAM111A (FAM111A)

<b>Product Code</b>	CSB-YP853484HU
<b>Abbreviation</b>	FAM111A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96PZ2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSCKKQRSRK HSVNEKCNMK IEHYFSPVSK EQQNNCSTSL MRMESRGDPR ATTNTQAQRF HSPKKNPEDQ TMPQNRTIYV TLKVNHRRNQ DMKLLKLTSE NSSLYMALNT LQAVRKEIET HQGQEMLVRG TEGIKEYINL GMPLSCFPEG GQVVITFSQS KSKQKEDNHI FGRQDKASTE CVKFYIHAIG IGKCKRRIVK CGKLHKKGRK LCVYAFKGET IKDALCKDGR FLSFLENDW KLIENNDTIL ESTQPVDELE GRYFQVEVEK RMVPSAAASQ NPSEKRNTC VLREQIVAQY PSLKRESEKI IENFKKKMKV KNGETLFELH RTTFGKVTKN SSSIKVVKLL VRLSDSVGYL FWDSATTGYA TCFVFKGLFI LTCRHVIDSI VGDGIEPSKW ATIIGQCVRV TFGYEELKDK ETNYFFVEPW FEIHNEELDY AVLKLENGQ QVPMELYNGI TPVPLSGLIH IIGHPYGEKK QIDACAVIPQ GQRAKQCQER VQSKKAESPE YVHMYTQRSF QKIVHNPDI TYDTEFFFGA SGSPVFDSKG SLVAMHAAGF AYTYQNETRS IIEFGSTMES ILLDIKQRHK PWYEEVFNQ QDVEMMSDED L
<b>Source</b>	Yeast
<b>Target Names</b>	FAM111A
<b>Protein Names</b>	Recommended name: Protein FAM111A
<b>Expression Region</b>	1-611
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.