



Recombinant Human DEP domain-containing protein 7 (DEPDC7)

Product Code	CSB-BP853487HU
Abbreviation	DEPDC7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QD5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MATVQEKAAA LNLSALHSPA HRPPGFSVAQ KPFGATYVWS SIINTLQTQV EVKKRRHRLK RHNDCEVVGSE AVDVIFSHLI QNKYFGDVDI PRAKVVRVCQ ALMDYKVFEA VPTKVFGKDK KPTFEDSSCS LYRFTTIPNQ DSQLGKENKL YSPARYADAL FKSSDIRSAS LEDLWENLSL KPANSPHVNI SATLSPQVIN EVWQEETIGR LLQLVDLPLL DSLLKQQEAV PKIPQPKRQS TMVNSSNYLD RGILKAYSDS QEDEWLSAAI DCLEYLPDQM VVEISRFPPE QPDRTDLVKE LLFDAIGRY YSSREPLLNHL SDVHNGIAEL LVNGKTEIAL EATQLLLKLL DFQNREEFRR LLYFMAVAAN PSEFKLQKES DNRMVVKRIF SKAIVDNKNL SKGKTDLLVL FLMDHQKDVF KIPGTLHKIV SVKLMAIQNG RDPNRDAGYI YCQRIDQRDY SNTEKTTKD ELLNLLKTL DEDSKLSAKEK KLLGQFYKC HPDIFIEHFG D
Source	Baculovirus
Target Names	DEPDC7
Protein Names	Recommended name: DEP domain-containing protein 7 Alternative name(s): Protein TR2/D15
Expression Region	1-511
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.