



Recombinant Human Rab-3A-interacting protein (RAB3IP)

Product Code	CSB-MP853488HU
Abbreviation	RAB3IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QF0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGLKMKGLS YDEAFAMAND PLEGFHEVNL ASPTSPDLLG VYESGTQEQT TSPSVIYRPH PSALSSVPIQ ANALDVSELP TQPVYSSPRR LNCAEISSIS FHVTDPAPCS TSGVTAGLTK LTTRKDNINA EREFLQGATI TEACDGSDDI FGLSTDLSLR LRSPSVLEVR EKGYERLKEE LAKAQRELKL KDEECERLSK VRDQLGQELE ELTASLFEEA HKMVREANIK QATAEKQLKE AQGKIDVLQA EVAALKTLVL SSSPTSPTQE PLPGGKTPFK KGHTRNKSTS SAMSGSHQDL SVIQPIVKDC KEADLSLYNE FRLWKDEPTM DRTCPFLDKI YQEDIFPCLT FSKSELASAV LEAVENNTLS IEPVGLQPIR FVKASAVECG GPKKCALTGQ SKSCKHRIKL GDSSNYYYIS PFCRYRITSV CNFFTYIRYI QQGLVKQQDV DQMFWEVMQL RKEMSLAKLG YFKEEL
Source	Mammalian cell
Target Names	RAB3IP
Protein Names	Recommended name: Rab-3A-interacting protein Short name= Rab3A-interacting protein Alternative name(s): Rabin-3 SSX2-interacting protein
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.