



Recombinant Human Ras-related protein Rab-40C (RAB40C)

Product Code	CSB-MP853500HU
Abbreviation	RAB40C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96S21
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGSQGSPVKS YDYLLKFLLV GDSDVGKGEI LESLQDGAAE SPYAYSNGID YKTTTILLDG RRVKLELWDT SGQGRFCTIF RSYSRGAQGI LLVYDITNRW SFDGIDRWIK EIDEHAPGVP RILVGNRLHL AFKRQVPTEQ ARAYAENKNCM TFFEVSPLCN FNVIESFTEL SRIVLMRHGM EKIWRPNRVF SLQDLCCRAI VSCTPVHLID KLPLPVTIKS HLKSFMSMANG MNAVMMHGRS YSLASGAGGG GSKGNSLKRS KSIRPPQSPP QNCSRSNCKI S
Source	Mammalian cell
Target Names	RAB40C
Protein Names	Recommended name: Ras-related protein Rab-40C Alternative name(s): Rar-like protein Ras-like protein family member 8C SOCS box-containing protein RAR3
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.