



# Recombinant Human Centrosomal protein of 89 kDa (CEP89)

<b>Product Code</b>	CSB-YP853505HU
<b>Abbreviation</b>	CEP89
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96ST8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLLGFRRGRR SHFKHIIHGL LPAASVAPKA AVPRTPPPRS PNPSPERPRS ALAAAILATT LTGRTVAIPQ PRQRSRSESD VSSVEQDSFI EPYATTSQLR PRPNWQSEMG RRSSLPSFET LDYGDEEDIE TQLSSSGKEL GDVSAREDRG GHSDDLAYAV HRNQVPLLHE VNSED DENIS HQDGFPGSPP APQRTQQKDG KHPVLNLKDE KPPLCEKPPP SPDITGRARQ RYTEITREKF EALKEENMDL NNMNQSLTLE LNTMKQAMKE LQLKLGMEK EKRKLKEAEK ASSQEVAAPE LLYLRKQAQE LVDENDGLKM TVHRLNVELS RYQTKFRHLS KEESLNIEGL PSKGPIPPWL LDIKYLSP LL LAYEDMMKEK DELNATLKEE MRMFRMRVQE VVKENEELHQ ELNKSSAVTS EEWRQLQTQA KLVLEENKLL LEQLEIQQRK AKDSHQERLQ EVSKLTKQLM LLEAKTHGQE KELAENREQL EILRAKQCEL KTHSDGKIAV EVHKSIVNEL KSQLQKEEEK ERAEMEELME KLTVLQAQKK SLLLEKNSLT EQNKALEAEL ERAQKINRKS QKKIEVLKKQ VEKAMGNEMS AHQYLANLVG LAENITQERD SLMCLAKCLE SEKDGVLNKV IKSNIIRLGKL EEKVKGYKKQ AALKLGDISH RLLEQQEDFA GKTAQYRQEM RHLHQVLKDK QEVLDQALQQ NREMEGELEV IWESTFREN RIRELLQDTL TRTGVDNPR ALVAPSLNGV SQADLLDGCD VCSYDLKSHA PTC
<b>Source</b>	Yeast
<b>Target Names</b>	CEP89
<b>Protein Names</b>	Recommended name: Centrosomal protein of 89 kDa Alternative name(s): Coiled-coil domain-containing protein 123
<b>Expression Region</b>	1-783
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.