



# Recombinant Schizosaccharomyces pombe Uncharacterized WD repeat-containing protein C32H8.09 (pi015, SPBC32H8.09)

<b>Product Code</b>	CSB-YP853522SXV
<b>Abbreviation</b>	pi015, SPBC32H8.09
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96WW0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDFTALFASL NPTFASVSHC GNWIASLSRS GHVLIRNSET LELHHVFLLN AQFIQKVYYL LWKPNLGAEK CHQICVASVD KVFVLDIVQH DYYASIQCDQ DPLSSISWSP SGELLWSSF DSKITVWSLN TQKGYLLPHV KTNVSKVYAL HPSMQFCTIL SRFNGSDCLQ FYQISKKAWI LLKECKLPTI DSTGIHWSPD GNWLAVLENV LDAVVYIYHR TGLLFHEYRP NRLIEVGFSD FEWSPFGKYL TLCYHDSTL HLETKTFSI VFRLHHCLQY TNTDLEMHIW EEKETIYEQQ MTYQKVHKLK TDFPEPSFCS ASKIRFNCDE TYAATITSKY PNVLWLWNLQ NKKLHTVLIQ KHHIVYFEWH PGRPDLVVIQ TKIRKESKIP SNATFLYFWA LSWNTPRVVG VPKKGFNIQK VQWLQPSEFS SRPVVICGE DAYTVAYIVE DEDESFT EVT QTIIISQEVLD NELETTQTIS IPS
<b>Source</b>	Yeast
<b>Target Names</b>	pi015
<b>Protein Names</b>	Recommended name: Uncharacterized WD repeat-containing protein C32H8.09
<b>Expression Region</b>	1-483
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.