



Recombinant Rhizobium loti UPF0169 lipoprotein Mll1543 (mll1543)

Product Code	CSB-MP853892RCU
Abbreviation	mll1543
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q98KC1
Product Type	Recombinant Protein
Immunogen Species	Rhizobium loti (strain MAFF303099) (Mesorhizobium loti)
Purity	>85% (SDS-PAGE)
Sequence	CMSSEKDID LSKYVDQTEP ADVLYNQGLA NLNAGRLDEA SKKFDAVDRQ HPYSEWARKS MVMGAFADYR KGSYDEAISS AKRYLALYPS TDDAPYAQYI IGLSYYRQIK DVTQDQKEAR QTLQTMQDLV TRWPTSEYVD DAKEKIRFAN DQLAGKEMQI GRYYLERREY IAAVKRFRTV VENYSNTRHV EEALARLTES YYAMGLTSEA QTAAAVLGTN YPDSPWYKDS YKLLQSNGLS PRENAGSWIS KAGKLITGA
Source	Mammalian cell
Target Names	bamD
Protein Names	Recommended name: UPF0169 lipoprotein Mll1543
Expression Region	32-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.