



Recombinant Human Protein FAM43A (FAM43A)

Product Code	CSB-MP854028HU
Abbreviation	FAM43A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N2R8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLPWKKHKFE LLAEAPPRQA SKPKGYAVSL HYSALSSLAR ACPEGALSRV GSMFRSKRKK LHITSEDPTY TVLYLGNATT IQARGDGCTD LAVGKIWSKS EAGRQGTKMK LTVSAQGIRM VHAERALRR PGHLYLLHRV TYCVADARLP KVFAWVYRHE LKHKAVMLRC HAVLVSKPEK AQAMALLLYQ TSANALAEFK RLKRRDDARH QQQELVGAHT IPLVPLRKLL LHGPCCYKPP VERSRSAPKL GSITEDLLGE QLEQELQEEE EEEQPEGCPE EEENRAAEGD PAEEEEAEQR ALVVAMHFEC GDLLDTLENG RGEALGGGGG SLGPGAGPPP LLLGSASDMK AELSQLISDL GELSGNDVTR TLQADLRVTR LLSGDSTGSE SSIEGGGPDA TSATAGDSSR QADGASADEP HSG
Source	Mammalian cell
Target Names	FAM43A
Protein Names	Recommended name: Protein FAM43A
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.