



Recombinant Human RING1 and YY1-binding protein (RYBP)

Product Code	CSB-MP854036HU
Abbreviation	RYBP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N488
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTMGDKKSPT RPKRQAKPAA DEGFWDCSV TFRNSAEAFK CSICDVRKGT STRKPRINSQ LVAQQVAQQY ATPPPKKEK KEKVEKQDKE KPEKDKEISP SVTKKNTNKK TKPKSDILKD PPSEANSIQS ANATTKTSET NHTSRPRLKN VDRSTAQQLA VTVGNVTVII TDFKEKTRSS STSSSTVTSS AGSEQQNQSS SGSESTDKGS SRSSTPKGDM SAVNDESF
Source	Mammalian cell
Target Names	RYBP
Protein Names	Recommended name: RING1 and YY1-binding protein Alternative name(s): Apoptin-associating protein 1 Short name= APAP-1 Death effector domain-associated factor Short name= DED-associated factor YY1 and E4TF1-associated factor 1
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.