



Recombinant Human Zinc-binding alcohol dehydrogenase domain-containing protein 2 (ZADH2)

Product Code	CSB-YP854039HU
Abbreviation	ZADH2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N4Q0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLRLVPTGAR AIVDMSYARH FLDFQGSAIP QAMQKLVVTR LSPNFREAVT LSRDCPVPLP GDGDLVLRNR FVGVNASDIN YSAGRYDPSV KPPFDIGFEG IGEVVALGLS ASARYTVGQA VAYMAPGSFA EYTVVPASIA TPVPSVKPEY LTLVSGTTA YISLKELGGL SEGKKVLVTA AAGGTGQFAM QLSKKAKCHV IGTCSSEDEKS AFLKSLGCDR PINYKTEPVG TVLKQEYPEG VDVVYESVGG AMFDLAVDAL ATKGRLIVIG FISGYQTPTG LSPVKAGTLP AKLLKKSASV QGFFLNHYLS KYQAAMSHLL EMCVSGDLVC EVDLGDLSPE GRFTGLESIF RAVNYMYMGK NTGKIVVELP HSVNSKL
Source	Yeast
Target Names	ZADH2
Protein Names	Recommended name: Zinc-binding alcohol dehydrogenase domain-containing protein 2 EC= 1.-.-
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.