



Recombinant Human Protein enabled homolog (ENAH)

Product Code	CSB-MP854072HU
Abbreviation	ENAH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N8S7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSEQSICQAR AAVMVYDDAN KKWVPAGGST GFSRVHIYHH TGNNTFRVVG RKIQDHQVVI NCAIPKGLKY NQATQTFHQW RDARQVYGLN FGSKEDANVF ASAMMHALEV LNSQETGPTL PRQNSQLPAQ VQNGPSQEEL EIQRRLQLEQ QRQKELERER LERERMERER LERERLERER LERERLEQEQ LERERQERER QERLERQERL ERQERLERQE RLDREERQERQ ERERLERLER ERQERERQEQ LEREQLEWER ERRISSAAP ASVETPLNSV LGDSSASEPG LQAASQPAET PSQQGIVLGP LAPPPPPPLP PGPAQASVAL PPPPGPPPPP PLPSTGPPPP PPPPPLPNQV PPPPPPPAP PLPASGFFLA SMSEDNRPLT GLAAAIAGAK LRKVSRMEDT SFPSGGNAIG VNSASSKTDG GRGNGPLPLG GSGLMEEMSA LLARRRRIAE KGSTIETEQA EDKGEDSEPV TSKASSTSTP EPTRKPWERT NTMNGSKSPV ISRRDSPRKN QIVFDNRSYD SLHRPKSTPL SQPSANGVQT EGLDYDRLKQ DILDREMRKEL TKLKEELIDA IRQELSKSNT A
Source	Mammalian cell
Target Names	ENAH
Protein Names	Recommended name: Protein enabled homolog
Expression Region	1-591
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.