



Recombinant Human Coiled-coil domain-containing protein 75 (CCDC75)

Product Code	CSB-MP854074HU
Abbreviation	CCDC75
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N954
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEEEDYMSD SFINVQEDIR PGLPMLRQIR EARRKEEKQQ EANLKNRQKS LKEEEQERRD IGLKNALGCE NKG FALLQKM GYKSGQALGK SGGGIVEPIP LNIKTGKSGI GHEASLKRKA EEKLESYRKK IHMKNQAE EK AAEQFRMRLK NKQDEM KLEG DLRRSQRACQ QLDVQKNIQV PREAWYWLRL EEETEEDEEE KEQDEDEYKS EDLSVLEKLQ ILTSYLREEH LYCIWCGTAY EDKEDLSSNC PGPTSADHD
Source	Mammalian cell
Target Names	GPATCH11
Protein Names	Recommended name: Coiled-coil domain-containing protein 75 Alternative name(s): G patch domain-containing protein 11
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.