



# Recombinant Human Uncharacterized protein C1orf127 (C1orf127)

<b>Product Code</b>	CSB-EP854077HU-B
<b>Abbreviation</b>	C1orf127
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N9H9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKCPMLRSRL GQESVHCGPM FIQVSRPLPL WRDNRQTPWL LSLRGELVAS LEDASLMGLY VDMNATTVTV QSPRQGLLQR WEVSGGQQAL PGVSFQPESE VLVHIPKQRL GLVKRGSYIE ETLSLRFLRV HQSNIFMVTE NKDFVVVSIP AAGVLQVQRC QEVGGTPTGTQ AFYRVDLSLE FAEMAAPVLW TVESFFQCVG SGTESPASTA ALRTPSPPS PGPETPPAGV PPAASSQVWA AGPAAQEWLS RDLLHRPSDA LAKKGLGPFL QTAKPARRGQ TSASILPRVV QAQRGPQPPP GEAGIPGHPT PPATLPSEPV EGVQASPWRP RPVLPHPAL TLPVSSDASS PSPPAPRPER PESLLVSGPS VTLTEGLGTV RPEQDPAKSP GSPLLLRLGLS SGDVAAPEPI MGEPGQASEE FQPLARPWRA TLAEEELVSH RSPGEPQETC SGTEVERPRQ TGPGLPREGA RGHMDLSSSE PSQDIEGPGL SILPARDATF STPSVRQPDP SAWLSSGPEL TGMPRVRLAA PLAVLPMEPL PPEPVRPAAL LTPEASSVGG PDQARYLESA PGWPVGQEEW GVAHTSSPPS TQTLSLWAPT GVLLPSLVEL EYPFQAGRGA SLQQELTEPT LALSAESHRP PELQDSVEGL SERPSR
<b>Source</b>	E.coli
<b>Target Names</b>	C1orf127
<b>Protein Names</b>	Recommended name: Uncharacterized protein C1orf127
<b>Expression Region</b>	1-656
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.