



Recombinant Human Zinc finger and SCAN domain-containing protein 4 (ZSCAN4)

Product Code	CSB-BP854087HU
Abbreviation	ZSCAN4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NAM6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MALDLRTIFQ CEPSENNLGS ENSAFQQSQG PAVQREEGIS EFSRMVLNSF QDSNNSYARQ ELQRLYRIFH SWLQPEKHSK DEIISLLVLE QFMIGGHEND KASVKEKWKS SGKNLERFIE DLTTDSINPP ALVHVHMQQG EALFSEDMPL RDVIVHLTKQ VNAQTTREAN MGTPSQTSQD TSLETGQGYE DEQDGNWSSS KTTRVNENIT NQGNQIVSLI IIQEENGRPR EEGGVSSDNP YNSKRAELVT ARSQEGSING ITFQGVPMVM GAGCISQPEQ SSPESALTHQ SNEGNTCEV HQKGSBGVQK SYKCECPKV FKYLCHLLAH QRRHRNERPF VCPECQKGF QISDLRVHQI IHTGKKPFTC SMCKKSFESHK TNLRSHERIH TGEKPYTCPF CKTSYRQSST YHRHMRTHEK ITPSPVPSTP EAS
Source	Baculovirus
Target Names	ZSCAN4
Protein Names	Recommended name: Zinc finger and SCAN domain-containing protein 4 Alternative name(s): Zinc finger protein 494
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.