



Recombinant Human Homeodomain-interacting protein kinase 4 (HIPK4)

Product Code	CSB-BP854106HU
Abbreviation	HIPK4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NE63
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSTIQSETDC YDIIIEVLGKG TFGEVAKGWR RSTGEMVAIK ILKNDAYRNR IINKNELKLLH CMRGLDPEEA HVIRFLEFFH DALKFYLVFE LLEQNLFEFQ KENNFAPLPA RHIRTVTLQV LTALARLKEK AIIHADLKPE NIMLVDQTRC PFRVKVIDFG SASIFSEVRY VKEPYIQSRF YRAPEILLGL PFCEKVDVWS LGCVMAELHL GWPLYPGNNE YDQVRYICET QGLPKPHLLH AACKAHFFFK RNPHPDAANP WQLKSSADYL AETKVRPLER RKYMLKSLDQ IETVNGGSVA SRLTFPDREA LAEHADLKSM VELIKRMLTW ESHERISPSA ALRHPFVSMQ QLRSAHETTH YYQLSLRSYR LSLQVEGKPP TPVVAEDGT PYYCLAEKE AAGMGSVAGS SPFFREEKAP GMQRAIDQLD DLSLQEAGHG LWGETCTNAV SDMMVPLKAA ITGHHVPDSG PEPILAFYSS RLAGRHKARK PPAGSKSDSN FSNLIRLSQV SPEDDRPCRG SSWEEGEHLG ASAEPLAILQ RDEDGPNIDN MTMEAERPDP ELFDPSSECPG EWLSEPDCTL ESVRGPRAQG LPPRRSHQHG PPRGATSFLQ HVTGHH
Source	Baculovirus
Target Names	HIPK4
Protein Names	Recommended name: Homeodomain-interacting protein kinase 4 EC= 2.7.11.1
Expression Region	1-616
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.