



Recombinant Human E3 ubiquitin-protein ligase ZSWIM2 (ZSWIM2)

Product Code	CSB-MP854109HU
Abbreviation	ZSWIM2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEG5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLRRGYKASE RRRHLSEKLS WHQDQALSSS IYLLREMGPT GFLLREEEPE YMDFRVFLGN PHVCNCSTFP KGGELCKHIC WLLKFKLP RNHESALQLG LGEREISDLL RGIHRVQTPQ PGTNDENEHV EEDGYIKQKE IDSEDICSIC QELLLEKKLP VTFCRFGCGN SIHIKCMKIL ANYQSTSNTS MLKCPLCRKE FAPKLILEE FKNSSKLVAE AEKERLDKHL GIPCNNCKQF PIEGKCYKCT ECIEYHLCQE CFDSCCHLSH TFTFREKRNRQ KWRSLEKRAD EVVKYIDTKN EIEEKMSHFQ EKQGQVYTPK HIVRSLPLQL ITKNSKLLAP GYQCCLLCLKA FHLGQHTRLL PCTHKFHRKC IDNWLFFHKCN SCPIDGQVIY NPLTWKNSAV NGQAHQSVSN RDIHLSKQK EPDLFIPGTG LVLKQNRLLGI LPSIPQCNFD ELNTPQSPKD AYENTTIDNL CSIKLDNSNS KKLTYDYKIS QHFPRYLQDL PTVSFGKIPS QTLLPPIVHK NIVCPTAMES PCISGKFHTS LSRMTKGCKC NNHNLKTPA TKIREDNKRS TLLPEDFNLI VNWSTAKLSL SKRYSNCMGE ITRKCSHLSR QPVSHSVNTK STELSLIEG VQL
Source	Mammalian cell
Target Names	ZSWIM2
Protein Names	Recommended name: E3 ubiquitin-protein ligase ZSWIM2 EC= 6.3.2.- Alternative name(s): MEK1-related protein X Short name= MEX ZZ-type zinc finger-containing protein 2 Zinc finger SWIM domain-containing protein 2
Expression Region	1-633
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.